Cyclonic Coriolis® technology qualification according to ISO 14698-1

Context

Bertin Technologies has developed a range of air samplers based on a patented innovative technology: Coriolis® cyclonic technology. In order to validate their physical and biological efficiency, these equipments have been tested in controlled rooms. Thus the technological performances of the whole range of Coriolis® products has been justified. The comparisons and protocols have been made using techniques described in ISO 14698-1 (Cleanrooms and associated controlled environments – biocontamination control).

Material

- Coriolis®MS, flow rate 600L/min
- Coriolis®μ, flow rate 300L/min
- Suspension of *Bacillus atrophaeus* spores
- Nebuliser
- Slit samplers
- Petri dishes (TS, TSA)
- Coriolis® collection liquid


Protocol

- Physical efficiency theoretical determination depending on the equipment design characteristics (air flow rate / cone diameter)
- Sampling in controlled environmental room: HPA (UK), WIS (GE), CEB (FR)
- Analysis by microbiological culture
- Comparison of the collection outputs with the reference method (slit samplers)

Results: Theoretical and experimental collection efficiency

1. Comparison with the theory (dry physical efficiency)
2. Experimental performances:


- Coriolis®μ – WIS: 69 % for 2 μm (9 samples)
- Coriolis®MS - HPA : 49 % for 0,8 μm (10 samples/diameter)
84 % for 4,6 μm
92 % for 10 μm
- Coriolis®MS - CEB : 50 % for 0,9 μm (3 samples/diameter)
57 % for 1,5 μm

Conclusions

- Experimental results (HPA, WIS, CEB) in accordance with the efficiency theory
- Equipments based on a single technology, all with an equal bioaerosols collection efficiency
- Coriolis® technology validated according to ISO 14698-1 and appropriate for cleanrooms and associated controlled environments applications

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