FEPAS® SAMPLE PREPARATION INSTRUCTIONS

Version 11 - Page 1 of 2, 28 September 2011

Please use the relevant instructions for individual samples supplied TREAT THE CONTENTS OF EACH CONTAINER AS A WHOLE – do not sub-sample

Instructions marked with * provide you with your 1/10 (-1) dilution

10 gram equivalent For the enumeration of target organism in beef/chicken/salad. *	Add 10 mL Buffered Peptone Water to rehydrate the sample. Shake the sample and stand at room temperature for 30 minutes.
	The sample can now be treated as if it were a routine 10 gram sample.
	Ensure that you rinse sample container with part of 90mL diluent before adding remaining diluent.
<u>10 gram</u>	Treat the sample as if it were a routine 10 gram sample.
For the enumeration of target organism in milk powder. *	Empty contents of vial into a Stomacher bag.
	Rinse sample container with part of 90mL diluent before adding remaining diluent.
	Stand for 30 minutes at room temperature.
<u>10 gram</u>	Treat the sample as if it were a routine 10 gram sample.
For the enumeration of target organism in animal feed/rice/flour.	However, be sure to use all the contents of the container and rinse the container thoroughly.
10 gram infant formula	Treat the sample as if it were a routine 10 gram sample.
For the detection of <i>Cronobacter sakazakii</i>	Ensure you treat the contents of the container as a whole –do not sub-sample. Rinse the container thoroughly.
25 gram equivalent	Add 20 mL Buffered Peptone Water to rehydrate the sample.
For the detection of target organism in chicken/beef/salad/fish.	Shake the sample and stand at room temperature for 30 minutes.
	The sample can now be treated as if it were a routine 25 gram sample.
	Rinse sample container with part of 225mL diluent before adding remaining diluent.
25 gram equivalent	Add 25 mL Buffered Peptone Water to rehydrate the sample.
For the detection of target organism in fresh milk/liquid egg * NB. these samples are supplied dried.	Shake the sample and stand at room temperature for 30 minutes.
	The sample can now be treated as if it were a routine 25 gram sample.
	Rinse sample container with part of 225mL diluent before adding remaining diluent.
25 gram For the detection of target organism in milk powder. *	Treat the sample as if it were a routine 25 gram sample.
	Empty all of the contents of vial into a Stomacher bag.
	 Rinse sample container with part of 225mL diluent before adding remaining diluent.
	Stand for 30 minutes at room temperature.

FEPAS® SAMPLE PREPARATION INSTRUCTIONS

Version 10 - Page 2 of 2, 28 September 2011

Please use the relevant instructions for individual samples supplied TREAT THE CONTENTS OF EACH CONTAINER AS A WHOLE – do not sub-sample

Instructions marked with * provide you with your 1/10 (-1) dilution

25 gram

For the **detection** of target organism in animal feed/chocolate/ cocoa powder/ground pepper/soft cheese.

- Treat each sample as if it were a routine 25 gram sample.
- Rinse sample container with part of 225mL diluent before adding remaining diluent.